## Sennheiser Hearing

## **TV** Hearing

Designed to give delightful hearing experience and effortless understanding when listening to media (TV, phone, tablet, computer, car)



Bluetooth (10m) / Airstream (15m)

Bluetooth (10m) / Airstream (15m)

**AUDIOLOGY & EAR CARE** 

Improve hearing in various situations where it

can be difficult to distinguish speech over

Bluetooth (10m) / Airstream (15m)

**Conversation Hearing** 

\*Plus shipping where applicable

Digital Wireless (70m)

Wireless Range

**Price\*** 

Weight

## Sennheiser Hearing



## Additional information

Ambient Awareness	Use Ambient Awareness to hear what's happening around you while you use your Sennheiser earbuds.
Tap Control	<ul> <li>Multi-function button on the earbud for convenient control</li> </ul>
Passive Noise Cancellation	<ul> <li>A secure in-ear design &amp; fit minimises external noise for distraction free listening</li> </ul>
Active Noise Cancelling	<ul> <li>Blocks out noise to make it easier to relax and focus</li> </ul>
<i>Automatic Scene Detection</i>	<ul> <li>Let speech enhancement match your surroundings with automatic scene detection</li> <li>Analyses the surrounding and automatically adjusts the Noise Canceller, directional microphones (beamforming) and volume to provide maximum enhancement when in Communication or Streaming mode</li> <li>In Communication mode the Classifier constantly monitors the environment and automatically adjusts the settings for the best speech understanding <ul> <li>Quiet: Noise canceller low; microphone beamformer off (= most natural speech)</li> <li>Moderate: Noise canceller medium; fixed microphone beamformer (= greatest sensitivity to the front)</li> <li>Loud: Noise canceller on; adaptive microphone beamformer (= greatest sensitivity to the front and additional cancellation to the loudest noise source)</li> </ul> </li> </ul>
Clarity modes	<ul> <li>Up to 20 dB high frequency gain for streamed signals</li> </ul>
<i>Conversation Clear</i> app	<ul> <li>For use with the ConC 400</li> <li>Simple guided profile setup using embedded sound samples to set preferred volume and clarity</li> <li>Tailor the performance of your earbuds in the Relax, Communicate, and Stream modes <ul> <li>Communication mode: see environmental noise level (soft, medium, loud); adjust speech volume; enable Clarity Boost</li> <li>Relax mode: general quite mode</li> <li>Streaming mode: balance transition between streaming and ambient volume around you</li> </ul> </li> </ul>
TV Clear app	<ul> <li>For use with the TV Clear and TVS 200</li> <li>Embedded setup guides</li> <li>Speech clarity selection, including demo song</li> <li>Volume control of ambient noise and streamed signal</li> <li>Configuration of controls, spoken messages etc</li> </ul>